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Brochure

This Brochure is additional and is not intended as a substitute for the Guarantee Booklet.

Manufacturer: B&T S.p.A. via Due Ponti 9, Forlì – Italy

Intended use of the product: Anti-decubitus mattress intended to optimise pressures during rest in order to reduce the risks of prolonged capillary occlusion and decubitus ulcer in individuals with reduced mobility.

Dorelan pressure-reducing mattresses in composite foam (hereinafter referred to simply as mattresses) are suitable to optimise pressures during rest in order to reduce the risks of prolonged capillary occlusion in individuals with reduced mobility. This characteristic is particularly appreciated in the health sector since it plays an important role in preventing pressure ulcers of stage one, if combined with specific changes of position. The reduction of pressure in the support zones is obtained by avoiding a steady posture which excessively stimulates sensitive body zones, if necessary by using accessories (e.g.: pillows).

The product complies with the essential safety requirements and with the UE Directive 93/42/CEE on medical devices, which was incorporated into Italian decree no. 46/97 of 24/02/97 and integration thereto, including the directive 2007/47/CE.

The medical device complies with UNI CEI EN ISO 14971 "Medical devices – Application of risk management to medical device" and UNI EN 12182 "Assistive products for person with disability – General requirement and test method". It also complies with requirements and specific tests according to the product type.

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INSTRUCTIONS ON THE CORRECT USE OF THE PRODUCT

It is recommended that the bed base slats are inserted into the frame by means of suitable holders, which allow the slats to be flexible and adjustable, thus providing flexibility to the whole base, including the perimeter.

It is recommended to check the adequacy of the product (i.e. mattress height) by placing a hand between the mattress and the slatted base, in order to check if the body touches the base (unsuitable support) or if it doesn't (suitable support).

Protect the mattress with disposable drawsheets which are air-permeable but impermeable to fluids.

The mattress is not fire-resistant.

The mattress is not sterile and cannot be sterilized.

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B&T S.p.A.
Tura Diano

B&T S.p.A. declares under sole responsibility that the products described in the table are in compliance with **Directive 93/42/EEC**.

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The registration on the website of Department of Health of the products listed above is available at the following address:

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